



CERTIFICATE OF ANALYSIS

Continuum™ Transfection Reagent (Catalog Number 400-700)

Lot Number: F122041

Storage Temperature: 2-8°C

Product Expiry: Oct2023

For Cell Culture, further Manufacturing or Research use only. Not for Direct Therapeutic use.

| Test | Methodology | Specification | Analysis |
|--------------------------------|-------------|---------------|-----------|
| Microbiological Testing | | | |
| Sterility | | | |
| Bacteria | | No Growth | No Growth |
| Fungi | | No Growth | No Growth |
| Biological Testing | | | |
| Functional Assay* | | ≥ 80% | Pass |

*Note: The functional assay consists of transfection of NIH3T3 cells with GFP expressing vector (pEGFP-N3). Flow cytometry analysis 24 hours post transfection must show at least 80% of GFP positive cells.

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Sofia Tann
Name

QC Lead
Title

04NOV2022
Date