



CERTIFICATE OF ANALYSIS
Recombinant Human/Mouse TGF- β 3
Transforming Growth Factor-Beta 3
Liquid
(Catalog Number 300-606)

Lot Number: G117045

Date Received: _____

Storage Temperature: 4°C

Product Expiry*: _____

For Cell culture, Further Manufacturing or Research use Only. Not for Direct Therapeutic Use.

*Expiration Date is one year from the date of receipt when stored at 4°C

Test	Methodology	Analysis
Biological Testing Endotoxin	LAL	0.069 EUs/ μ g
Biological Testing Purity	Reducing and Non-Reducing SDS PAGE	>95%
Biological Activity ED ₅₀	Inhibition of IL-4 induced HT-2cell proliferation	5.5 x 10 ⁶ units/mg

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Alina Z Khan
Name

QC Associate I
Title

17 MAR 2022
Date