



## CERTIFICATE OF ANALYSIS

### **Yssel's Custom T-cell Medium Without Phenol Red (Catalog Number 900-133)**

Lot Number: M122083

Date of Manufacture: Apr2022

Storage Temperature:  $\leq -10^{\circ}\text{C}$

Product Expiry: Apr2024

For Cell culture, Further Manufacturing or Research use Only. Not for Direct Therapeutic Use.

**Product description:** Yssel's Custom T-cell Medium without Phenol Red is sterile-filtered through a 0.2  $\mu\text{m}$  filter prior to freeze.

Test	Methodology	Specification	Analysis
<b>Biological Testing</b>			
Endotoxin	USP<85>, EP 2.6.14	<5 EU/mL	<0.500 EU/mL
<b>Microbiological Testing</b>			
Sterility	USP<71>, EP 2.6.1		
Bacteria		No Growth	No Growth
Fungi		No Growth	No Growth
<b>Physical Testing</b>			
Osmolality	USP<785>, EP 2.2.35	250-350 mOsm/Kg	297 mOsm/kg
pH	USP<791>	6.5-7.5	7.0

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Julie Wilson  
Name

QA Specialist  
Title

13 MAY 2022  
Date