



CERTIFICATE OF ANALYSIS

Yssel's Custom T-cell Medium Without Phenol Red (Catalog Number 900-133)

Lot Number: M122087

Date of Manufacture: Apr2022

Storage Temperature: $\leq -10^{\circ}\text{C}$

Product Expiry: Apr2024

For Cell culture, Further Manufacturing or Research use Only. Not for Direct Therapeutic Use.

Product description: Yssel's Custom T-cell Medium without Phenol Red is sterile-filtered through a 0.2 μm filter prior to freeze.

| Test | Methodology | Specification | Analysis |
|--------------------------------|---------------------|-----------------|-------------|
| Biological Testing | | | |
| Endotoxin | USP<85>, EP 2.6.14 | <5 EU/mL | <0.05 EU/mL |
| Microbiological Testing | | | |
| Sterility | USP<71>, EP 2.6.1 | | |
| Bacteria | | No Growth | No Growth |
| Fungi | | No Growth | No Growth |
| Physical Testing | | | |
| Osmolality | USP<785>, EP 2.2.35 | 250-350 mOsm/Kg | 294 mOsm/Kg |
| pH | USP<791> | 6.5-7.5 | 7.1 |

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

VANESSA YOUNG

Name

QA ASSOCIATE II

Title

12 MAY 2022

Date