


Product:

Bovine Serum Albumin (BSA) Solution—30%
 Ultra-High Monomer Reagent Grade, Salt Adjusted

Catalog #: 700-114
Lot #: C54J00L
Origin: United States (100%)
Manufacture Date: Apr 2019
Expiration Date: Apr 2022

TEST	METHODOLOGY	SPECIFICATION	RESULTS
Microbiological Testing			
Microbial Count	USP <61>	< 1 cfu/mL	< 1 cfu/mL
Physical Testing			
pH	Neat	6.8 -7.3	7.1
Biochemical Testing			
Protein	Chemistry Analyzer	28 to 31%	30%
Protease	Enzymatic Assay	< 0.005 Units/mL	< 0.005 Units/mL
Sodium Chloride	Chloride Analyzer	0.6 to 0.7 g/dL	0.6 g/dL
Bovine igG	Chemistry Analyzer	None Detected	None Detected
NEFA	Enzymatic Assay	FIO	<0.001%

This material was prepared from bovine plasma obtained from animals that have received ante- and post-mortem health inspections from APHIS Inspectors within facilities inspected by the USDA.


 29 Jul 2020

Created by/Date


 29 JUL 2020

Verified by/Date