



Product: **Continuum™ Transfection Reagent**

Catalog #: 400-700
Lot #: F00T72J

Manufacture Date: Dec 2018
Expiration Date: Dec 2020

TEST	SPECIFICATION	RESULTS
Microbiological Testing		
Sterility		
Bacteria	Not Detected	Not Detected
Fungi	Not Detected	Not Detected
Biological Testing		
Functional Assay*	≥80%	Pass

*Note: The functional assay consists of transfection of NIH3T3 cells with GFP vector (pEGFP-N3) Flow Cytometry analysis 24 hours post transfection must show at least 80% of GFP positive cells.

For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.

 09 Apr 2019

Created by / Date



Verified by / Date