



**Product:** **Continuum™ Transfection Reagent**

**Catalog #: 400-700**  
**Lot #: F03T72H**

**Manufacture Date: Jul 2016**  
**Expiration Date: Oct 2019**

TEST	SPECIFICATION	RESULTS
<b>Microbiological Testing</b>		
Sterility		
Bacteria	Not Detected	Not Detected
Fungi	Not Detected	Not Detected
<b>Biological Testing</b>		
Functional Assay*	≥80%	Pass

\*Note: The functional assay consists of transfection of NIH3T3 cells with GFP vector (pEGFP-N3) Flow Cytometry analysis 24 hours post transfection must show at least 80% of GFP positive cells.

*For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.*

Erin Jones 11 Apr 2019  
Created by / Date

M. M. M. M. 12 Apr 2019  
Verified by / Date