




Product: Continuum™ Transfection Reagent

Catalog #: 400-700
Lot #: F08S72H

Manufacture Date: Jul 2016
Expiration Date: Apr 2019

TEST	SPECIFICATION	RESULTS
Microbiological Testing		
Sterility		
Bacteria	Not Detected	Not Detected
Fungi	Not Detected	Not Detected
Biological Testing		
Functional Assay*	≥80%	Pass

*Note: The functional assay consists of transfection of NIH3T3 cells with GFP vector (pEGFP-N3) Flow Cytometry analysis 24 hours post transfection must show at least 80% of GFP positive cells.



Created by



Verified by