



**Product:** Continuum™ Transfection Reagent

**Catalog #: 400-700**

**Lot #: F121073**

**Expiration Date: Nov 2022**

TEST	SPECIFICATION	RESULTS
<b>Microbiological Testing</b>		
Sterility		
Bacteria	Not Detected	Not Detected
Fungi	Not Detected	Not Detected
<b>Biological Testing</b>		
Functional Assay*	≥80%	Pass

\*Note: The functional assay consists of transfection of NIH3T3 cells with GFP vector (pEGFP-N3) Flow Cytometry analysis 24 hours post transfection must show at least 80% of GFP positive cells.

*For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.*

 11 NOV 2021  
 \_\_\_\_\_  
 Created by / Date

 11 NOV 2021  
 \_\_\_\_\_  
 Verified by / Date