



Product: **Human Fibronectin**
0.1% Solution (Stabilized)
0.2µm sterile-filtered

Catalog #: 400-135
Lot #: F49R45H

Manufacture Date: Oct 2016
Expiration Date: Oct 2018

<i>TEST</i>	<i>METHODOLOGY</i>	<i>SPECIFICATION</i>	<i>RESULTS</i>
Biological Testing			
Endotoxin	LAL	Report	<4 ng/mg
Viral Testing			
HBsAg		Non-Reactive	Non-Reactive
Anti-HIV-1/HIV-2		Non-Reactive	Non-Reactive
Syphilis		Negative	Negative
HIV-NAT		Non-Reactive	Non-Reactive
HCV-NAT		Non-Reactive	Non-Reactive
Physical Testing			
Protein		Report	4 mg
Purity	SDS Electrophoresis	Report	>90%



 Created by


 Verified by