



**GeminiBio**  
**CERTIFICATE OF ANALYSIS**

**GemCell™ Select Human Serum AB**  
**(Part Number 100-812)**

**Lot Number:**

H121015

**Date of Manufacture:** MAR 2021

**Product Expiry:** MAR 2026

**Origin:** United States

**Storage Temperature:** ≤ -10°C

For Cell culture, Further Manufacturing or Research use Only. Not for Direct Therapeutic Use.

**Product description:** *GemCell™ Select Human Serum AB is collected from healthy male donors of the AB serotype at FDA-licensed facilities located in the United States. Donor units are tested for infectious disease markers prior to processing and found to be non-reactive. Gemcell™ Select Human Serum is converted to serum from human plasma using bovine thrombin and sterile-filtered through a 0.1 µm filter prior to freeze.*

Test	Methodology	Specification	Analysis
<b>Biological Testing</b>			
Endotoxin	USP<85>, EP 2.6.14	< 10.0 EU/ mL	<0.750 EU/mL
Hemoglobin	Fleming, AF and Woolf, AJ (1985)	<20.0 mg/ dL	6.4 mg/dL
<b>Microbiological Testing</b>			
Sterility	USP<71>, EP 2.6.1		
Bacteria		No Growth	No Growth
Fungi		No Growth	No Growth
Mycoplasma	Barile, MF and Kern, J (1971), EP 2.6.7	Not Detected	Not Detected
<b>Viral Testing</b>			
HBsAg	ABBOTT ChLIA	Non-Reactive	Non-Reactive
Anti-HCV	ABBOTT ChLIA	Non-Reactive	Non-Reactive
Anti-HIV-1/ HIV-2	ABBOTT ChLIA	Non-Reactive	Non-Reactive
Syphilis	ASI RPR	Negative	Negative
HBV-NAT	ROCHE NAT/ PCR	Not Detected	Not Detected
HIV-NAT	ROCHE NAT/ PCR	Not Detected	Not Detected
HCV-NAT	ROCHE NAT/ PCR	Not Detected	Not Detected
HAV RNA		Non-Reactive	Non-Reactive
Parvo B19 DNA		Non-Reactive	Non-Reactive
Atypical Antibodies		Negative	Negative
Viral Testing Finished Pool Level			
HBcAB		Non-Reactive	Non-Reactive
HTLV-I/II		Non-Reactive	Non-Reactive
West Nile Virus		Not Detected	Not Detected
Chagas		Non-Reactive	Non-Reactive
HEV IgG		Non-Reactive	Non-Reactive



**GeminiBio**  
**CERTIFICATE OF ANALYSIS**

**GemCell™ Select Human Serum AB**  
**(Part Number 100-812)**

Test	Lot Number: Methodology	H121015 Specification	Analysis
<b>Physical Testing</b>			
Osmolality	USP<785>, EP 2.2.35	260 – 350 mOsm/ kg	311 mOsm/ kg
pH @ RT	USP<791>	Test and Report	8.14
<b>Biochemistry</b>			
Albumin		Test and Report	3.4 g/dL
ALT (SGPT)		Test and Report	9 U/L
AST (SGOT)		Test and Report	12 U/L
Bilirubin, Total		Test and Report	0.2 mg/dL
BUN		Test and Report	13 mg/dL
Calcium		Test and Report	>15.0 mg/dL
Chloride		Test and Report	110 mmol/L
Cholesterol		Test and Report	128 mg/dL
Creatinine		Test and Report	0.72 mg/dL
Glucose		Test and Report	92 mg/dL
Phosphorus		Test and Report	2.9 mg/dL
Potassium		Test and Report	3.7 mg/dL
Protein, Total	Biuret	Test and Report	5.1 g/dL
Sodium		Test and Report	>165 mmol/L
Triglycerides		Test and Report	99 mg/dL
Uric Acid		Test and Report	3.1 mg/dL

*GemCell Select™ Human Serum AB is collected from a maximum of 16 healthy male donors of the AB serotype at FDA-licensed facilities in the United States. This material is defibrinated from source and the viral testing is performed on the individual donor units. All other testing is performed on the final product pool prior to release. The bovine thrombin that is used as part of the conversion to serum is sourced from controlled herds located in the United States and the source cattle are ante and post-mortem inspected by a U.S. Veterinary Service Inspector where they were deemed free of infectious and contagious diseases. All animals used in the production of the thrombin were from a natural beef program in accordance with FDA regulations. Material is derived from human blood and should be considered biohazardous. Universal precautions should be used when handling this material. Precipitate may develop in this product upon freezing and/or thawing; this occurrence does not impact culture performance.*

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Henna Z Khan  
Name

QC Associate I  
Title

24 MAY 2021  
Date