



## CERTIFICATE OF ANALYSIS

### Addendum

### Human Serum AB

### (Off the Clot)

### (Catalog Number 100-318)

Lot Number: H122043

Date of Manufacture: Dec2022

Product Expiry: Dec2027

Origin: United States

Storage Temperature:  $\leq -10^{\circ}\text{C}$

For Cell culture, Further Manufacturing or Research use Only. Not for Direct Therapeutic Use.

**Product description:** Human serum AB is collected from healthy male donors of the AB serotype at FDA-licensed facilities located in the United States. Donor units are tested for infectious disease markers prior to processing and found to be non-reactive. Human Serum AB(off the clot) is sterile-filtered through a 0.1  $\mu\text{m}$  filter prior to freeze.

| Test   | Methodology                    | Analysis    |
|--|--------------------------------|-------------|
| Viral Testing Finished Pool<br>Donor HTLV I/II Antibody Screen |                                | Nonreactive |
| <u>Franklin Quiroz</u><br>Name                                 | <u>Karen Zuniga</u><br>Name    |             |
| <u>QC lead</u><br>Title  | <u>QA Associate I</u><br>Title |             |
| <u>01 MAR 2023</u><br>Date                                     | <u>01 MAR 2023</u><br>Date     |             |