



CERTIFICATE OF ANALYSIS
Human Recombinant Insulin
10mg/mL Solution
(Catalog Number 800-112)

Lot Number: H423008

Date of Manufacture: Apr2023

Storage Temperature: ≤-10°C

Product Expiry: Apr2026

For Cell Culture, further Manufacturing or Research use only. Not for Direct Therapeutic use.

Test	Methodology	Specification	Analysis
Microbiological Testing			
Sterility	USP<71>		
Bacteria		No Growth	No Growth
Fungi		No Growth	No Growth
Physical Testing			
pH @ RT	USP<791>	Report	2.95

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Karen Zuniga

Name

QA Associate I

Title

01MAY2023

Date