



CERTIFICATE OF ANALYSIS
Phosphate - Buffered Saline Solution
(Catalog Number 600-108)

Lot Number: M122044

Date of Manufacture: Mar2022

Product Expiry: Mar2024

Storage Temperature: Ambient

Product description: *Phosphate - Buffered Saline Solution is sterile-filtered through a 0.2 µm filter prior to storage.*

Test	Methodology	Specification	Analysis
Biological Testing			
Endotoxin	USP <85>	Report	<0.050 EU/mL
Microbiological Testing			
Sterility	USP <71>, EP 2.6.1		
Bacteria		No Growth	No Growth
Fungi		No Growth	No Growth
Physical Testing			
pH@RT	USP<791>	Report	7.45
Osmolality	USP<785>	Report	283 mOsm/kg

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Karen Zuniga

QA Associate I

Name 28 MAR 2022

Title

Date