



**CERTIFICATE OF ANALYSIS**  
**Phosphate - Buffered Saline Solution**  
**(Catalog Number 600-108)**

Parent Lot Number: M63S00L

Date of Manufacture: May2020

Sub-Lot Number: M63S00L-18, M63S00L-19, M63S00L-20, M63S00L-21, M63S00L-22, M63S00L-23, M63S00L-24, M63S00L-25, M63S00L-26, M63S00L-27

Product Expiry: May2022

For Transport of Specimens only.  
Not to be taken internally.

For In Vitro Diagnostic Use  
Storage Temperature: Store at 15 - 30°C

**Product description:** Phosphate - Buffered Saline Solution is sterile-filtered through a 0.2 µm filter prior to storage.

Test	Methodology	Specification	Analysis
<b>Biological Testing</b>			
Endotoxin	USP <85>	Report	<0.500 EU/mL
<b>Microbiological Testing</b>			
37°C Sterility Check		No Growth	No Growth
<b>Physical Testing</b>			
pH@RT	USP<791>	Report	7.38
Osmolality	USP<785>	Report	284 mOsm/Kg

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Megha Khatali

Name

QA Associate III

Title

17 JUN 2020

Date