



CERTIFICATE OF ANALYSIS
Phosphate - Buffered Saline Solution
(Catalog Number 600-108)

Parent Lot Number: M72S00L

Date of Manufacture: Jun2020

Sub-Lot Number: M72S00L-1, M72S00L-2, M72S00L-4,
M72S00L-5, M72S00L-6

Product Expiry: Jun2022

For Transport of Specimens only.
Not to be taken internally.

For In Vitro Diagnostic Use
Storage Temperature: Store at 15 - 30°C

Product description: *Phosphate - Buffered Saline Solution is sterile-filtered through a 0.2 µm filter prior to storage.*

Test	Methodology	Specification	Analysis
Biological Testing			
Endotoxin	USP <85>	Report	<0.500
Microbiological Testing			
37°C Sterility Check		No Growth	No Growth
Physical Testing			
pH@RT	USP<791>	Report	7.44
Osmolality	USP<785>	Report	286

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Chunty Francisco
Name
25 SEP 2020
Date

QC Associate I
Title