



## CERTIFICATE OF ANALYSIS

### Saline Solution

(Catalog Number 600-107)

Parent Lot Number: M76S00L

Date of Manufacture: Jul2020

Sub-Lot Number: M76S00L-1, M76S00L-2, M76S00L-3, M76S00L-4, M76S00L-5, M76S00L-6, M76S00L-7, M76S00L-8, M76S00L-9, M76S00L-10, M76S00L-11, M76S00L-12, M76S00L-13, M76S00L-14, M76S00L-15, M76S00L-16, M76S00L-17, M76S00L-18, M76S00L-19, M76S00L-20, M76S00L-21

Product Expiry: Jul2022

For Transport of Specimens only.  
Not to be taken internally.

For In Vitro Diagnostic Use  
Storage Temperature: Store at 15 - 30°C

Product description: Saline Solution is sterile-filtered through a 0.2 µm filter prior to storage.

Test	Methodology	Specification	Analysis
<b>Biological Testing</b>			
Endotoxin	USP <85>	Report	<0.500 EU/mL
<b>Microbiological Testing</b>			
37°C Sterility Check		No Growth	No Growth
<b>Physical Testing</b>			
pH@RT	USP<791>	Report	6.17
Osmolality	USP<785>	Report	268 mOsm/Kg

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Megha Khawasi  
Name

QA Associate III  
Title

11 Aug 2020  
Date