



CERTIFICATE OF ANALYSIS
Saline Solution
(Catalog Number 600-107)

Parent Lot Number: M80S00L

Date of Manufacture: Aug2020

Sub-Lot Number: M80S00L-1, M80S00L-2, M80S00L-3, M80S00L-4, M80S00L-5, M80S00L-6, M80S00L-7, M80S00L-8, M80S00L-9, M80S00L-10, M80S00L-11, M80S00L-12, M80S00L-13, M80S00L-14, M80S00L-15, M80S00L-16, M80S00L-17, M80S00L-18, M80S00L-19, M80S00L-20, M80S00L-21, M80S00L-22, M80S00L-23, M80S00L-24, M80S00L-25

Product Expiry: Aug2022

For Transport of Specimens only.
Not to be taken internally.

For In Vitro Diagnostic Use
Storage Temperature: Store at 15 - 30°C

Product description: Saline Solution is sterile-filtered through a 0.2 µm filter prior to storage.

Test	Methodology	Specification	Analysis
Biological Testing			
Endotoxin	USP <85>	Report	<0.625 EU/mL
Microbiological Testing			
37°C Sterility Check		No Growth	No Growth
Physical Testing			
pH@RT	USP<791>	Report	5.65
Osmolality	USP<785>	Report	269 mOsm/Kg

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Megha Khairani
Name
26 Aug 2020
Date

QA Associate III
Title