



CERTIFICATE OF ANALYSIS

Saline Solution

(Catalog Number 600-107)

Parent Lot Number: M85S00L

Date of Manufacture: Aug2020

Sub-Lot Number: M85S00L-1, M85S00L-2, M85S00L-3, M85S00L-4, M85S00L-5, M85S00L-6, M85S00L-7, M85S00L-8, M85S00L-9, M85S00L-10, M85S00L-11, M85S00L-12, M85S00L-13, M85S00L-14, M85S00L-15, M85S00L-16, M85S00L-17, M85S00L-18, M85S00L-19, M85S00L-20, M85S00L-21, M85S00L-22, M85S00L-23, M85S00L-24, M85S00L-25

Product Expiry: Aug2022

For Transport of Specimens only.
Not to be taken internally.

For In Vitro Diagnostic Use
Storage Temperature: Store at 15 - 30°C

Product description: Saline Solution is sterile-filtered through a 0.2 µm filter prior to storage.

Test	Methodology	Specification	Analysis
Biological Testing			
Endotoxin	USP <85>	Report	<0.500 EU/mL
Microbiological Testing			
37°C Sterility Check		No Growth	No Growth
Physical Testing			
pH@RT	USP<791>	Report	5.25
Osmolality	USP<785>	Report	280 mOsm/Kg

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Megha Khalasi
Name
10 Sep 2020
Date

QA Associate III
Title