



CERTIFICATE OF ANALYSIS
Saline Solution
(Catalog Number 600-107)

Parent Lot Number: M86S00L

Date of Manufacture: Aug2020

Sub-Lot Number: M86S00L-1, M86S00L-2, M86S00L-3, M86S00L-4, M86S00L-5, M86S00L-6

Product Expiry: Aug2022

For Transport of Specimens only.
Not to be taken internally.

For In Vitro Diagnostic Use
Storage Temperature: Store at 15 - 30°C

Product description: *Saline Solution is sterile-filtered through a 0.2 µm filter prior to storage.*

Test	Methodology	Specification	Analysis
Biological Testing			
Endotoxin	USP <85>	Report	<0.591 EU/mL
Microbiological Testing			
37°C Sterility Check		No Growth	No Growth
Physical Testing			
pH@RT	USP<791>	Report	6.05
Osmolality	USP<785>	Report	268 mOsm/Kg

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Megha Khatali

Name

10 Sep 2020

Date

QA Associate III

Title