



## CERTIFICATE OF ANALYSIS

### Saline Solution

(Catalog Number 600-107)

Parent Lot Number: M89S00L

Date of Manufacture: Sep2020

Sub-Lot Number: M89S00L-1, M89S00L-2, M89S00L-3, M89S00L-4, M89S00L-5, M89S00L-6, M89S00L-7, M89S00L-8, M89S00L-9, M89S00L-10, M89S00L-11

Product Expiry: Sep2022

For Transport of Specimens only.  
Not to be taken internally.

For In Vitro Diagnostic Use  
Storage Temperature: Store at 15 - 30°C

Product description: Saline Solution is sterile-filtered through a 0.2 µm filter prior to storage.

Test	Methodology	Specification	Analysis
<b>Biological Testing</b>			
Endotoxin	USP <85>	Report	<0.500 EU/mL
<b>Microbiological Testing</b>			
37°C Sterility Check		No Growth	No Growth
<b>Physical Testing</b>			
pH@RT	USP<791>	Report	6.09
Osmolality	USP<785>	Report	271 mOsm/kg

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Megha Khalasi

Name

24 SEP 2020

Date

QA Associate III

Title