



## CERTIFICATE OF ANALYSIS

### Addendum rev2

#### **GemPure Select Water for Injection (WFI) Quality Water (Catalog Number 901-001-021, 901-001-101, 901-001-201)**

Lot Number: M122106

Date of Manufacture: Jun2022

Expiration Date: Jun2023

Storage Temperature: 2°C – 30°C

For Cell culture, Further Manufacturing or Research use Only. Not for Direct Therapeutic Use.

**Product Description:** *GemPure Select Water for Injection (WFI) Quality Water meets the USP monograph for Water for Injection packaged in bulk for commercial use.*

Test	Methodology	Specification	Analysis
Appearance	Visual Inspection	Clear and colorless liquid	Clear and colorless liquid
<b>Biological Testing</b>			
Endotoxin	USP <85>, EP 2.6.14	≤ 0.25 EU/mL	<0.01 EU/mL
<b>Microbiological Testing</b>			
Sterility	USP <71>, EP 2.6.1		
Bacteria		No Growth	No Growth
Fungi		No Growth	No Growth
<b>Biochemistry</b>			
Total Organic Carbon (Before Packaging)	USP <643>	≤ 500 ppb	7.84 ppb
Total Organic Carbon (Packaged)	USP <643>	≤ 8000 ppb	4410 ppb
<b>Physical Testing</b>			
pH @ 25°C	USP <791>	5.0 – 7.0	5.9
Conductivity @ 25°C (Packaged)	USP <645>	≤ 5 μS/cm	1.276 μS/cm
Conductivity @ 25°C (Before Packaging)	USP <645>	≤ 1.3 μS/cm	0.09 μS/cm

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and it has been confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Jequeline Wilson  
Name

QA Specialist  
Title

16 MAY 2024  
Date