



CERTIFICATE OF ANALYSIS

Human Albumin 25% Solution

(Catalog Number 800-120)

Lot Number: H424003

Date of Manufacture: Jan2024

Product Expiry: May2030

Origin: United States

Storage Temperature: $\leq 10^{\circ}\text{C}$

Test	Methodology	Specification	Analysis
Microbiological Testing			
Sterility	USP<71>		
Bacteria		No Growth	No Growth
Fungi		No Growth	No Growth
Viral Testing			
	21 CFR 610.40		
HBsAg		Non-Reactive	Non-Reactive
Anti-HCV		Non-Reactive	Non-Reactive
Anti-HIV-1/HIV-2		Non-Reactive	Non-Reactive
Syphilis		Negative	Negative
HIV-1 NAT		Non-Reactive	Non-Reactive
HBV-NAT		Non-Reactive	Non-Reactive
HCV-NAT		Non-Reactive	Non-Reactive
Physical Testing			
Purity		Test and Report	97%
Protein		Test and Report	25.0 g/dL
Endotoxin	USP<85>	Test and Report	0.523 EU/mL
pH	USP<791>	Test and Report	6.96

For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use. Results shown were obtained by carefully performed methods believed to be reliable. However, since some results may vary for specific tests depending upon methodology and other variables, it is suggested that tests for which results are particularly important be repeated by the user of this product.

All human plasma is collected from stringently screened donors at FDA-licensed collection centers in the United States.

VANESSA YOUNG
Name

QA SPECIALIST
Title

03JUL2024
Date