



## CERTIFICATE OF ANALYSIS

### Addendum

### Human Serum AB

### Off-the-Clot

(Catalog Number 100-318-100, 100-318-005)

Lot Number: H424011

Date of Manufacture: Apr2024

Product Expiry: Apr2029

Origin: United States

Storage Temperature:  $\leq -10^{\circ}\text{C}$

For Cell Culture, further Manufacturing or Research use only. Not for direct Therapeutic use.

**Product Description:** Human Serum AB Off-the-Clot is collected from healthy male donors of the AB serotype at FDA-licensed facilities located in the United States. Donor units are tested for infectious disease markers prior to processing and found to be non-reactive. Human Serum AB Off-the-Clot is sterile-filtered through a  $0.1\ \mu\text{m}$  filter prior to freezing.

Test	Methodology	Analysis
<b>Viral Testing Finished Pool</b>		
Donor HTLV I/II Antibody Screen		Nonreactive
<u>Korrina Nicewonger</u>	<u>VANESSA YOUNG</u>	
Name	Name	
<u>QA Associate II</u>	<u>QA SPECIALIST</u>	
Title	Title	
<u>13 JUN 2024</u>	<u>13 JUN 2024</u>	
Date	Date	